

# PLANETARY DISC MILL

Code : M434



- Used whenever highest fineness is required
- Robust model with 4 grinding stations
- Short grinding time assured by centrifugal force of planetary mixing
- Jars are made of hardened steel (also available in stainless steel upon request)
- Used for grinding minerals, ores, coal, coke, alloys, bentonite, bones, cement clinkers, concrete, glass, gypsum, limestone, polymers, seeds ... etc.
- Feed size is < 10 mm
- Final fineness is < 10 microns
- Speed control for required applications
- Automatic cut-off energy if the cover is opened

## **TECHNICAL SPECIFICATIONS**

- Robust model with 4 grading stations
- Feed size is < 10 mm
- Final fineness is < 10 microns
- Speed control for required applications
- Automatic cut-off energy if the cover is opened
- 220 - 240 V / 50 Hz

**ORDERING INFORMATION**

Item	Code
PLANETARY DISC MILL	M434X04XU
GRINDING CELL for PLANETARY DISC MILL	M434P001H