

BALL MILL

Code : M438



- Used to dry-grind lime stones, ceramic, glass and slag

TECHNICAL SPECIFICATIONS

- Drum made of stainless steel (300 x 300 mm)
- Adjustable Rotation Speed (Max. 70 rpm)
- Supplied with:
- Grinding Balls:
 - 43 x Ø38.1 mm
 - 67 x Ø31.75 mm
 - 10 x Ø25.4 mm
 - 71 x Ø19.05 mm
 - 94 x Ø15.87 mm

- Grinding Mills:
 - 10 x 20 mm
 - 10 x 25 mm
 - 10 x 30 mm
- 220-240 V/ 50-60 Hz

ORDERING INFORMATION

Item	Code
BALL MILL	M438X00XU
ABRASIVE CHARGES SET for M438	M438P001U