

# SKID RESISTANCE TESTER

Code : R190



- Used for measuring surface frictional properties. It is used to measure the energy loss when a rubber slider edge is propelled over a test surface.
- Suitable for laboratory as well as field tests on flat surfaces, and for polish value measurements on curved laboratory specimens from accelerated polishing-wheel tests.

## **STANDARDS**

ASTM E-103 • BS 812-144 • EN 1097

## **SUPPLIED WITH**

- 6 rubber sliders
- Specimen Holder

**ORDERING INFORMATION**

Item	Code
SKID RESISTANCE TESTER	R190X00XH