

FALLING HEAD PERMEABILITY SET

Code : T250



- Used to determine the permeability of fine-grained soils

STANDARDS

EN 17892-11

SUPPLIED WITH

- Permeameter stand with 3 tubes
- Plated-Steel Falling-head permeability cell - Ø 100 mm (1 ea)
- Water Reservoir Tank

ORDERING INFORMATION

Item	Code
FALLING HEAD PERMEABILITY SET	T250X00XH
MANOMETER STAND with 4 TUBES (Falling Head)	T250P001H
RESERVOIR TANK (Falling Head)	T250P002H
PERMEABILITY CELL (Falling Head) - Ø 4"	T250P003H
PERMEABILITY CELL (Falling Head) - Ø 6"	T250P004H