

UTM, ELECTROMECHANICAL

Code : U002



- Specially designed for testing components and for routine materials such as sheets, wires, strips, ropes, belts, screws, composites, nylon ribbon, rubber, geosynthetics, geotextiles ... etc
- Features precision guidance for moving crosshead and a low wearball screw trouble-free tensile tests.

STANDARDS

EN ISO 6892 • EN ISO 4136 • ASTM A370 • EN ISO 527-1

ORDERING INFORMATION

Item	Code
UTM, ELECTROMECHANICAL - 10 kN	U002A001U

UTM, ELECTROMECHANICAL - 50 kN	U002A005U
UTM, ELECTROMECHANICAL - 100 kN	U002A010U
UTM, ELECTROMECHANICAL - 150 kN	U002A015U
UTM, ELECTROMECHANICAL - 200 kN	U002A020U
UTM, ELECTROMECHANICAL - 300 kN	U002A030U
UTM, ELECTROMECHANICAL - 600 kN	U002A060U
UTM, ELECTROMECHANICAL - 750 kN	U002A075U
UTM, ELECTROMECHANICAL - 1000 kN	U002A100U
MECHANICAL GRIPPERS SET for U002xxxxx	U002P001H
BENDING APPARATUS for U002xxxxx	U002P002H

OTHER PHOTOS

