

VOLUMETRIC FLASK

Code : W003



- A volumetric flask in engineering testing laboratories is used for accurately preparing and diluting solutions to a specific volume. Its precision makes it ideal for measuring exact concentrations in chemical, environmental, or material testing, ensuring consistency and accuracy in experiments.

ORDERING INFORMATION

| Item | Code |
|--------------------------|-----------|
| VOLUMETRIC FLASK - 5 ml | W003G005H |
| VOLUMETRIC FLASK - 10 ml | W003G010H |
| VOLUMETRIC FLASK - 25 ml | W003G025H |
| VOLUMETRIC FLASK - 50 ml | W003G050H |

| | |
|----------------------------|-----------|
| VOLUMETRIC FLASK - 100 ml | W003G100H |
| VOLUMETRIC FLASK - 150 ml | W003G150H |
| VOLUMETRIC FLASK - 250 ml | W003G250H |
| VOLUMETRIC FLASK - 500 ml | W003G500H |
| VOLUMETRIC FLASK - 1000 ml | W003G01KH |
| VOLUMETRIC FLASK - 2000 ml | W003G02KH |
| VOLUMETRIC FLASK - 3000 ml | W003G03KH |
| VOLUMETRIC FLASK - 4000 ml | W003G04KH |
| VOLUMETRIC FLASK - 5000 ml | W003G05KH |
| VOLUMETRIC FLASK - 6000 ml | W003G06KH |

OTHER PHOTOS

