

CRUCIBLE TONG

Code : W022



- Crucible tongs in engineering testing laboratories are used to handle hot crucibles or containers, especially when they are heated to high temperatures. They provide a safe way to move and manipulate hot materials during processes like melting, burning, or other high-temperature chemical reactions.

ORDERING INFORMATION

| Item | Code |
|------------------------|-----------|
| CRUCIBLE TONG - 250 mm | W022M250H |
| CRUCIBLE TONG - 300 mm | W022M300H |
| CRUCIBLE TONG - 450 mm | W022M450H |