

## PORCELAIN MORTAR and PESTLE

Code : W026



- A porcelain mortar and pestle in engineering testing laboratories are used for grinding, crushing, and mixing solids into fine powders. The non-porous nature of porcelain ensures that no material is absorbed, making it ideal for preparing samples for chemical or material analysis without contamination.

### **ORDERING INFORMATION**

Item	Code
PORCELAIN MORTAR and PESTLE - Ø 80 mm	W026C080H
PORCELAIN MORTAR and PESTLE - Ø 100 mm	W026C100H
PORCELAIN MORTAR and PESTLE - Ø 130 mm	W026C130H
PORCELAIN MORTAR and PESTLE - Ø 160 mm	W026C160H

---

PORCELAIN MORTAR and PESTLE - Ø 190 mm	W026C190H
PORCELAIN MORTAR and RUBBER-HEADED PESTLE	W027P001H
GLASS MORTAR and PESTLE - Ø 100 mm	W028G100H
GLASS MORTAR and PESTLE - Ø 120 mm	W028G120H
GLASS MORTAR and PESTLE - Ø 150 mm	W028G150H
GLASS MORTAR and PESTLE - Ø 180 mm	W028G180H

**OTHER PHOTOS**

