

# CENTRIFUGE EXTRACTOR

Code : A031



- Used to determine the bitumen percentage in asphalt mixtures.
- The extractor is available in two capacities (1500 g and 3000 g), and is equipped with a speed control knob.
- *Filter Papers should be ordered separately*

## **STANDARDS**

ASTM D2172 • AASHTO T-164

## **TECHNICAL SPECIFICATIONS**

- Capacity:
  - 1500 g
  - 3000 g
- Speed control knob
- Electronic Brake for fast stopping
- 220 - 240 V / 50 - 60 Hz

### **ORDERING INFORMATION**

| Item                                         | Code      |
|----------------------------------------------|-----------|
| CENTRIFUGE EXTRACTOR - 1500 g                | A031M15XU |
| CENTRIFUGE EXTRACTOR - 3000 g                | A031M30XU |
| FILTER PAPER for A031M15XU                   | A031P001H |
| SPARE BOWL for CENTRIFUGE EXTRACTOR - 1500 g | A031P002H |
| FILTER PAPER for A031M30XU                   | A031P003H |
| SPARE BOWL for CENTRIFUGE EXTRACTOR - 3000 g | A031P004H |