

# CONCRETE PIPE TESTING MACHINE

Code : B008



- Supplied with a graphical LCD Data Acquisition Control System and an ability to store and report the results.
- The calibration for the equipment can be easily done on 5 pts.
- Pipe length and diameter should be specified at the time of inquiry.

## **STANDARDS**

ASTM C301 • ASTM C497 • EN 1916

## **TECHNICAL SPECIFICATIONS**

| Capacity (kN)      | 400        | 600        | 1000       |
|--------------------|------------|------------|------------|
| Pipe Diameter (mm) | 200 - 1500 | 200 - 2000 | 200 - 3700 |

|                  |            |            |             |
|------------------|------------|------------|-------------|
| Pipe Length (mm) | 500 - 1500 | 500 - 2000 | 1000 - 3000 |
|------------------|------------|------------|-------------|

## ORDERING INFORMATION

| Item                                    | Code      |
|-----------------------------------------|-----------|
| CONCRETE PIPE TESTING MACHINE - 400 kN  | B008A04XU |
| CONCRETE PIPE TESTING MACHINE - 600 kN  | B008A06XU |
| CONCRETE PIPE TESTING MACHINE - 1000 kN | B008A10XU |
| V-SHAPED SUPPORT - 1 m                  | B008P001H |