

CAPPING SET

Code : B015



- Comprises a melting pot and a capping apparatus.
- The melting pot is designed to melt capping compounds and features an adjustable temperature control via a thermoregulator.
- A built-in pilot lamp indicates the heating status.
- The pot is constructed with a double-walled frame for enhanced durability and insulation.

STANDARDS

ASTM C31 • ASTM C192 • EN 12390-3 • AASHTO T23

SUPPLIED WITH

- Capping apparatus (Used for Ø 100 mm and Ø 150 mm cylinders)

TECHNICAL SPECIFICATIONS

- Melting pot (4.5 lt):
 - Thermostatically Controlled
 - Insulated double-wall
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

| Item | Code |
|----------------------|-----------|
| CAPPING FRAME | B015P001H |
| MELTING POT | B015P002U |
| CYLINDER CAPPING SET | B015X00XU |

OTHER PHOTOS

