

# KELLY-BALL PENETRATION APPARATUS

Code : B022



- Used to determine the workability of concrete as an alternative to Slump Test, where 1" of kelly ball penetration corresponds to 2" of slump value
- Coated for rust protection
- Weight of the movable parts: 30 lb (13.6 kg)

## **STANDARDS**

ASTM C360

## **TECHNICAL SPECIFICATIONS**

- Metal hemispherical cylinder:
  - Diameter : Ø 6"

- Height : 4 5/8"

**SUPPLIED WITH**

- Metal hemispherical cylinder and graduated handle

**ORDERING INFORMATION**

| Item                 | Code         |
|----------------------|--------------|
| KELLY BALL APPARATUS | B 022 X00A H |