

VIBRATING TABLE

Code : B125



- Designed from robust steel and has a very steady frame
- Used to compact the concrete specimens in the laboratory
- Equipped with a vibrating motor
- Supplied with a clamp assembly for moulds.

STANDARDS

ASTM C192 • EN 12390-2

ORDERING INFORMATION

Item	Code
VIBRATING TABLE (Small)	B125X001U

VIBRATING TABLE (Large)	B125X002U
-------------------------	-----------

VIBRATING TABLE (Large) - 380V	B125X002V
--------------------------------	-----------

CLAMPING APPARATUS for VIBRATING TABLE (Small)	B125P001U
--	-----------

CLAMPING APPARATUS for VIBRATING TABLE (Large)	B125P002U
--	-----------