

# CEMENT TUBE SAMPLER

Code : C028



- The Cement Tube Sampler is specifically designed for obtaining representative samples of hydraulic cement from individual packages, ensuring accurate quality control assessments.
- Constructed from durable brass, the sampler maintains structural integrity and resists corrosion, providing consistent performance across numerous sampling operations.
- Its slender, cylindrical design allows for easy penetration into cement bags, facilitating the extraction of samples without significant disturbance to the remaining material.
- The sampler's dimensions are optimized to collect sufficient quantities of cement for comprehensive testing, aligning with international standards for material sampling.
- Equipped with a comfortable hardwood handle, the sampler offers ease of use during manual operations, reducing operator fatigue and enhancing efficiency.

## **STANDARDS**

ASTM C183 • EN 196-7 • AASHTO T127

**TECHNICAL SPECIFICATIONS**

- Material: Brass
- Outer Diameter: 32 mm
- Length: 700 mm
- Handle: Hardwood

**ORDERING INFORMATION**

Item	Code
CEMENT TUBE SAMPLER (Packaged)	C028X00XH
CEMENT SLOTTED TUBE SAMPLER (Bulk)	C029X00XH