

## **CRADLE**

Code : G100P001H



- Designed to suspend specimens beneath a weighing balance within a Specific Gravity Frame, facilitating accurate mass measurements in both air and water for specific gravity calculations.
- Constructed from high-quality stainless steel rod, ensuring durability and resistance to corrosion during frequent immersion processes.
- Compatible with ALFA's Specific Gravity Frame, seamlessly integrating into existing laboratory setups for aggregate testing.

### **STANDARDS**

ASTM C127 • AASHTO T85 • BS 812:2 • EN 12390-7 • EN 1097-6

### **ORDERING INFORMATION**

Item	Code
SPECIFIC GRAVITY FRAME	G100X00XH
CRADLE	G100P001H
DENSITY BASKET (Ø 200 x 200 mm / Aperture: 3.35 mm)	G100P002H
WATER TANK	G100P003H
DENSITY / SOUNDNESS BASKET (Ø 250 x 250 mm / Aperture: 1.00 mm)	G100P004H
DENSITY / SOUNDNESS BASKET (Ø 120 x 160 mm / Aperture: 3.35 mm)	G100P005H
SPECIFIC GRAVITY PYCNOMETER	G103X00XH