


VERNIER CALIPER

Code : G709



- Used to measure the dimensions of an object, generally by placing two movable points of the instrument across the object or span to be measured.
- Many types of calipers permit reading out a measurement on a ruled scale, a dial, or an electronic digital display. A common association is to calipers using a sliding vernier scale.

ORDERING INFORMATION

| Item | Code |
|--------------------------|-----------|
| DIGITAL CALIPER - 150 mm | G709X001B |
| DIGITAL CALIPER - 200 mm | G709X002B |
| DIGITAL CALIPER - 300 mm | G709X003B |

| | |
|----------------------------|-----------|
| DIGITAL CALIPER - 500 mm | G709X004B |
| VERNIER CALIPER - 150 mm | G709X005H |
| VERNIER CALIPER - 300 mm | G709X006H |
| VERNIER CALIPER - 500 mm | G709X007H |
| DIGITAL MICROMETER - 25 mm | G709X101B |

OTHER PHOTOS

