

MOISTURE TIN

Code : G720



- Moisture tins in engineering testing laboratories are used to hold soil or material samples during oven drying for moisture content determination. They are durable, heat-resistant, and help ensure accurate weight measurements before and after drying.

ORDERING INFORMATION

| Item | Code |
|---|-----------|
| MOISTURE TIN, with LID - Small - Ø 55 x 35 mm | G720X001H |
| MOISTURE TIN, with LID - Medium | G720X002H |
| MOISTURE TIN, with LID - Large | G720X003H |