

MAGNESIUM SULFATE

Code : G806



- Magnesium sulfate in engineering testing laboratories is often used in the durability testing of concrete and aggregates. It is used in the magnesium sulfate soundness test to assess the resistance of aggregates to weathering and cyclic environmental conditions, helping determine their suitability for construction.

ORDERING INFORMATION

Item	Code
TRICHLOROETHYLENE	G801X00XH
METHYLENE BLUE POWDER - 100 g	G802X00XH
SULPHUR	G803X00XH
GRAPHITE	G804X00XH
SODIUM SULFATE	G805X00XH

MAGNESIUM SULFATE	G806X00XH
BARIUM CHLORIDE	G807X00XH
CALCIUM CHLORIDE	G808X00XH
SODIUM HYDROXIDE	G809X00XH
FORMALDEHYDE (Stock Solution)	G810X00XH
SODIUM HEXAMETAPHOSPHATE	G811X00XH
SILICA GEL	G812X00XH
PARAFIN WAX	G813X00XH
KAOLINITE	G814X00XH
QUANTAB CHLORIDE TITRATOR 1175 TYPE	G815X00XH
QUANTAB CHLORIDE TITRATOR 1176 TYPE	G816X00XH