


SLAKE DURABILITY

Code : K130



- Used to determine the durability of the rock specimens to weakening and disintegrating when subjected to the climatic slaking effects.
- Equipped with a motorized drive unit that rotates two drums. Drums rotate at 20 rpm inside acrylic water tank.
- The drums are made of stainless steel mesh.

TECHNICAL SPECIFICATIONS

- Motorized Drive Unit
- Rotates 2 drums
- Rotates at 20 rpm
- Acrylic Water Tank
- Drums of Stainless Steel Mesh

- Drum Size : Ø 140 mm x 100 mm
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item	Code
SLAKE DURABILITY - 2 Drums	K130M02AU
SLAKE DURABILITY - 4 Drums	K130M04AU
DRUM for SLAKE DURABILITY	K130P001H
TANK for SLAKE DURABILITY	K130P002H