

SOIL > ATTERBERG LIMITS

SHRINKAGE LIMIT TEST SET

Code : T037



- Used to determine the maximum amount of moisture in the sample at which the soil will not shrink while drying

STANDARDS

ASTM D427 • AASHTO T92 • BS 1377:2 • NF P94-060-1

SET

Item

CRYSTALLIZING CUP (Ø 57 x 31 mm)

SHRINKAGE DISH (Ø 44 x 12 mm) [2 ea]

SHRINKAGE PRONG PLATE

EVAPORATING DISH [300 ml - Ø 125 mm]

MEASURING CYLINDER - 250 ml

SPATULA - 100 mm

CARRYING CASE for T037

ORDERING INFORMATION

Item	Code
SHRINKAGE LIMIT SET	T 037 X00X H
CRYSTALLIZING CUP (Ø 57 x 31 mm)	T 037 P001 H
SHRINKAGE DISH (Ø 44 x 12 mm)	T 037 P002 H
SHRINKAGE PRONG PLATE	T 037 P003 H
EVAPORATING DISH [300 ml - Ø 125 mm]	W 137 C300 H
MEASURING CYLINDER - 250 ml	W 052 G250 H
SPATULA - 100 mm	G 717 X100 H
CARRYING CASE for T037	T 037 P004 H