

## **SAMPLE EXTRUDER, Universal**

Code : U020



- Used to extract the sample from its mould. It can be used with proctor test moulds, CBR test moulds, Marshall test moulds and other moulds having the diameter 100 mm to 6".
- The machine is actuated using a manual hydraulic mechanism.

### **TECHNICAL SPECIFICATIONS**

- Manual-Hydraulic Type
- Used for Ø 100 mm (Ø 4") and Ø 150 mm (Ø 6") diameter samples
- Can be used with the moulds of Proctor, CBR, Marshall Tests

### **AVAILABLE MODELS**

| Code         | Item                                |
|--------------|-------------------------------------|
| U 020 H100 H | SAMPLE EXTRUDER, UNIVERSAL - 100 mm |
| U 020 H300 H | UNIVERSAL SAMPLE EXTRUDER - 300 mm  |
| U 020 H450 H | UNIVERSAL SAMPLE EXTRUDER - 400 mm  |
| U 020 H600 H | UNIVERSAL SAMPLE EXTRUDER - 600 mm  |
| U 020 M300 U | MOTORIZED SAMPLE EXTRUDER - 300 mm  |

### ORDERING INFORMATION

| Code         | Item                                |
|--------------|-------------------------------------|
| U 020 H100 H | SAMPLE EXTRUDER, UNIVERSAL - 100 mm |
| U 020 H300 H | UNIVERSAL SAMPLE EXTRUDER - 300 mm  |
| U 020 H450 H | UNIVERSAL SAMPLE EXTRUDER - 400 mm  |
| U 020 H600 H | UNIVERSAL SAMPLE EXTRUDER - 600 mm  |
| U 020 M300 U | MOTORIZED SAMPLE EXTRUDER - 300 mm  |
| U 020 P001 H | ADAPTOR for [U 020] - Ø 38 mm       |
| U 020 P002 H | ADAPTOR for [U 020] - Ø 50 mm       |
| U 020 P003 H | ADAPTOR for [U 020] - Ø 70 mm       |
| U 020 P004 H | ADAPTOR for [U 020] - Ø 106 mm      |

### OTHER PHOTOS

