

LABORATORY BOTTLE

Code : W006



- A laboratory bottle in engineering testing laboratories is used for storing, transporting, and handling liquids, chemicals, or samples. These bottles are typically made of glass or durable plastic and are designed to safely contain substances while providing ease of access and measurement during testing procedures.

ORDERING INFORMATION

Item	Code
LABORATORY BOTTLE - 25 ml	W006G025H
LABORATORY BOTTLE - 50 ml	W006G050H
LABORATORY BOTTLE - 100 ml	W006G100H

LABORATORY BOTTLE - 250 ml	W006G250H
----------------------------	-----------

LABORATORY BOTTLE - 500 ml	W006G500H
----------------------------	-----------

LABORATORY BOTTLE - 1000 ml	W006G01KH
-----------------------------	-----------

LABORATORY BOTTLE - 2000 ml	W006G02KH
-----------------------------	-----------

OTHER PHOTOS

