

# DESSICATOR

Code : W008



- A desiccator in engineering testing laboratories is used to store samples or materials in a dry environment by absorbing moisture from the air. It helps prevent sample contamination and ensures accurate results in tests that are sensitive to humidity, such as moisture content analysis or when working with hygroscopic materials.

## **ORDERING INFORMATION**

Item	Code
DESSICATOR - 200 mm	W008G296H
DESSICATOR - 250 mm	W008G344H
DESSICATOR - 300 mm	W008G420H
VACUUM DESSICATOR - 200 mm	W009G296H

---

VACUUM DESSICATOR - 250 mm

W009G344H

VACUUM DESSICATOR - 300 mm

W009G420H