

# VACUUM DESSICATOR

Code : W009



- A vacuum desiccator in engineering testing laboratories is used to store samples in a low-pressure, dry environment. By applying a vacuum, it removes moisture from the air and from the sample itself, ensuring that samples remain dry and uncontaminated during storage or preparation for sensitive tests.

## **ORDERING INFORMATION**

Item	Code
DESSICATOR - 200 mm	W008G296H
DESSICATOR - 250 mm	W008G344H
DESSICATOR - 300 mm	W008G420H
VACUUM DESSICATOR - 200 mm	W009G296H
VACUUM DESSICATOR - 250 mm	W009G344H

VACUUM DESSICATOR - 300 mm

W009G420H