

PYCNOMETER (Gay Lussac)

Code : W014



- A pycnometer (Gay Lussac type) in engineering testing laboratories is used to determine the specific gravity or density of liquids or powders. It consists of a precisely calibrated glass container, often with a stopper, ensuring accurate measurements by accounting for temperature variations during testing. This is essential for density and quality control in chemical, material, and environmental testing.

ORDERING INFORMATION

Item	Code
PYCNOMETER (Gay Lussac) - 5 ml	W014G005H
PYCNOMETER (Gay Lussac) - 10 ml	W014G010H
PYCNOMETER (Gay Lussac) - 25 ml	W014G025H

PYCNOMETER (Gay Lussac) - 50 ml	W014G050H
---------------------------------	-----------

PYCNOMETER (Gay Lussac) - 100 ml	W014G100H
----------------------------------	-----------

PYCNOMETER (Gay Lussac) - 250 ml	W014G250H
----------------------------------	-----------

PYCNOMETER (Gay Lussac) - 500 ml	W014G500H
----------------------------------	-----------

PYCNOMETER (Gay Lussac) - 1000 ml	W014G01KH
-----------------------------------	-----------

PYCNOMETER (Gay Lussac) - 2000 ml	W014G02KH
-----------------------------------	-----------

OTHER PHOTOS

