

## LE CHATELIER FLASK

Code : W019



- A Le Chatelier flask in engineering testing laboratories is used for studying equilibrium reactions and measuring the solubility of gases in liquids under different conditions. It is typically used in chemical analysis to observe the effect of temperature or pressure changes on the solubility and volume of gases, helping in studies of chemical equilibria.

### **ORDERING INFORMATION**

Item	Code
LE CHATELIER MOULD	C023X00EH
LE CHATELIER SOUNDNESS TEST SET	C025X00EH
LE CHATELIER BATH - 6 lt	G040X06XU
LE CHATELIER FLASK - 250 ml	W019G250H

---

LE CHATELIER STAND for [G 040 X06X U]

G040P001H

SPARE HEATER for LE CHATELIER BATH

G040P004H