

## PORCELAIN MORTAR and RUBBER-HEADED PESTLE

Code : W027



- A porcelain mortar and rubber-headed pestle in engineering testing laboratories are used for grinding and mixing materials without causing contamination or damage. The rubber-headed pestle is gentle on the mortar's surface and is ideal for handling delicate or fragile samples, making it useful for preparing samples in chemical, biological, or material testing.

### **ORDERING INFORMATION**

Item	Code
PORCELAIN MORTAR and PESTLE - Ø 80 mm	W026C080H
PORCELAIN MORTAR and PESTLE - Ø 100 mm	W026C100H
PORCELAIN MORTAR and PESTLE - Ø 130 mm	W026C130H
PORCELAIN MORTAR and PESTLE - Ø 160 mm	W026C160H

---

PORCELAIN MORTAR and PESTLE - Ø 190 mm	W026C190H
PORCELAIN MORTAR and RUBBER-HEADED PESTLE	W027P001H
GLASS MORTAR and PESTLE - Ø 100 mm	W028G100H
GLASS MORTAR and PESTLE - Ø 120 mm	W028G120H
GLASS MORTAR and PESTLE - Ø 150 mm	W028G150H
GLASS MORTAR and PESTLE - Ø 180 mm	W028G180H