

## GLASS MORTAR and PESTLE

Code : W028



- A glass mortar and pestle in engineering testing laboratories are used for grinding and mixing materials, particularly for substances that may react with porcelain or other materials. The smooth surface of glass ensures no contamination and is ideal for preparing fine powders or combining chemicals in precise quantities for testing.

### **ORDERING INFORMATION**

Item	Code
PORCELAIN MORTAR and PESTLE - Ø 80 mm	W026C080H
PORCELAIN MORTAR and PESTLE - Ø 100 mm	W026C100H
PORCELAIN MORTAR and PESTLE - Ø 130 mm	W026C130H

---

PORCELAIN MORTAR and PESTLE - Ø 160 mm	W026C160H
PORCELAIN MORTAR and PESTLE - Ø 190 mm	W026C190H
PORCELAIN MORTAR and RUBBER-HEADED PESTLE	W027P001H
GLASS MORTAR and PESTLE - Ø 100 mm	W028G100H
GLASS MORTAR and PESTLE - Ø 120 mm	W028G120H
GLASS MORTAR and PESTLE - Ø 150 mm	W028G150H
GLASS MORTAR and PESTLE - Ø 180 mm	W028G180H

**OTHER PHOTOS**

