

# PLASTIC MEASURING CYLINDER

Code : W032



- A plastic measuring cylinder in engineering testing laboratories is used for accurately measuring the volume of liquids. It is lightweight and durable, making it suitable for handling a variety of chemicals or materials, especially when glass containers might be too fragile or heavy. It provides precise measurements for experiments and sample preparation.

## **ORDERING INFORMATION**

Item	Code
PLASTIC MEASURING CYLINDER - 10 ml	W032P010H
PLASTIC MEASURING CYLINDER - 25 ml	W032P025H
PLASTIC MEASURING CYLINDER - 50 ml	W032P050H
PLASTIC MEASURING CYLINDER - 100 ml	W032P100H

---

PLASTIC MEASURING CYLINDER - 250 ml	W032P250H
PLASTIC MEASURING CYLINDER - 500 ml	W032P500H
PLASTIC MEASURING CYLINDER - 1000 ml	W032P01KH
PLASTIC MEASURING CYLINDER - 2000 ml	W032P02KH