

# PLASTIC VOLUMETRIC FLASK

Code : W033



- A plastic volumetric flask in engineering testing laboratories is used for accurately preparing and diluting solutions to a precise volume. It is lightweight and durable, making it suitable for handling chemicals and liquids, especially when glass flasks are prone to breakage. It ensures precision in chemical analysis and sample preparation.

## **ORDERING INFORMATION**

Item	Code
PLASTIC VOLUMETRIK FLASK - 25 ml	W033P025H
PLASTIC VOLUMETRIK FLASK - 50 ml	W033P050H
PLASTIC VOLUMETRIK FLASK - 100 ml	W033P100H
PLASTIC VOLUMETRIK FLASK - 250 ml	W033P250H

---

PLASTIC VOLUMETRIK FLASK - 500 ml

W033P500H

PLASTIC VOLUMETRIK FLASK - 1000 ml

W033P01KH